



Prestan[®] Ultralite[®] Manikins



The Ultralite Manikin is amazingly easy to set up and is sold in a convenient and lightweight 4-pack for efficient training-on-the-go. This durable manikin offers an affordable tool for CPR training, with all the quality and realism you expect from Prestan.

Prestan[®] Ultralite[®] Manikins



Features & Benefits:

- Highly compact, lightweight design – complete 4-pack weighs only 13 lb
- Unique clamshell opens easily for simple insertion of face shield lung bags, which attach to the manikin to prevent slipping
- Visible chest rise when ventilated
- Head tilt/chin lift teaches students to correctly open the airway
- One-piece, disposable face shield lung bag
- Realistic chest compressions to a full 2" depth
- Compatibility with AED trainer pads
- Life tested to a half-million compressions
- Latex-free
- Stackable torsos and heads fit in a compact carrying case and are quick and easy to assemble at training sites
- Molded, pliable skin is realistic to the eye and touch; resists dirt, grime, and grease; is durable and easy to clean; and allows AED pads to adhere to the manikin without leaving adhesive residue behind
- Realistic head features, including a prominent jawline for simulation of jaw thrust technique
- Realistic torso with prominent collarbones, an accurate nipple line, obvious xiphoid process, rib-cage lines, breastbone structure, and an upper chest sized for the easy application of AED pads



Ordering Information

- | | |
|---------|-------------------------------|
| #10-226 | Manikin Pack |
| #10-225 | Face shield lung bags (50/pk) |
| #10-224 | Face shields (50/pk) |

Includes

4 manikin torsos, 4 manikin heads, 4 compression pistons, nylon piston bag, 50-pack face-shield lung bags, nylon carrying case with detachable shoulder strap and instruction sheet.